



COLLAGEN RESEARCH CENTER

1-10-1 Kamikiyoto, Kiyose, Tokyo 204-0013, Japan
Tel +81-424-95-1995 Fax+81-424-95-1990

Product List

Type II Collagen

No.	Name	Quantity	MRP	Outline
K41	0.3% Type II Collagen (Liquid), Bovine	24mg/8ml	28,000	immunization
K42	Type II Collagen (Lyophilized powder), Bovine	24mg	26,000	immunization
K43	Type II Collagen (Lyophilized powder), Chicken	10mg	34,000	immunization
K44	Type II Collagen (Lyophilized powder), Porcine	10mg	34,000	immunization
K45	Type II Collagen (Liquid), Rat ELISA	200µg	34,000	ELISA
K46	Type II Collagen (Liquid), Bovine ELISA	100µg	24,000	ELISA
K47	Type II Collagen (Liquid), Chicken ELISA	100µg	24,000	ELISA
K48	Type II Collagen (Liquid), Porcine ELISA	100µg	24,000	ELISA
K49	Type II Collagen (Liquid), Human ELISA	100µg	34,000	ELISA
K49S	Type II Collagen (Lyophilized powder), Monkey ELISA	500µg	34,000	ELISA

Type IV Collagen

No.	Name	Quantity	MRP	Outline
K34	Human Placenta Type IV Collagen, lyophilized	1mg	34,000	The major band of molecular weight 70,000–80,000.
K35	Type IV Collagen NC1(Liquid), Bovine	1mg/ml	47,000	It can be used in the development and screening of renal disease therapeutics. It can also be used in autoimmune disease research. Induction of autoimmune Nephritis in Monkey (This is the first in the world.)

Type I Collagen

No.	Name	Quantity	MRP	Outline
K31	Type I Collagen	100mg/50ml	28,000	0.2% Type I collagen

Collagenase

No.	Name	Quantity	MRP	Outline
K10	Collagenase from tadpoles	10µl	28,000	Purified product derived from tadpoles.
K30	Collagenase, from bacteria	50ml	34,000	From aerobic bacterium (V. Hollisae). Cell culture, Study on improvement of tissue fibrosis.

Antibody

No.	Name	Quantity	MRP	Outline
K34POLY	Rabbit Anti Type IV collagen Antiserum	50 µl	47,000	Rabbit polyclonal to Collagen IV. This is the immune staining of tissue basement membrane
K35POLY	Rabbit Anti Type IV Collagen NC1 domain Antiserum	50 µl	47,000	
K35MONO	Anti Type IV Collagen NC1 domain Monoclonal Antibody	50 µl	47,000	K35MONO stain only failure renal, does not stain healthy renal and recovery renal in indirect immunofluorescence staining. (It was confirmed in monkey and human) Patent. US, EU, JP
K51	Anti Type II collagen antibody	500µl	28,000	"Specimen of anti-type II collagen antibody Antigen / bovine type II collagen Immunized animals / Lewis rats"
K52	Anti Type II collagen monoclonal antibody	20 µl	28,000	Antigen / chicken. Crossover / Bovine, Human, Porcine, Rat

KIT

No.	Name	Quantity	MRP	Outline
K11	Collageno Kit CLN100	100 tests	48,000	Collagenase activity assay kit with FITC-type I collagen substrate. Research for Rheumatoid
K12	Collageno Kit Type II	48 tests	57,000	Collagenase activity assay kit with FITC-type II collagen substrate. It can be used in research for bone metabolism and cartilage cells and joint disease.
K35KIT	Type IV Collagen NC1 domain Test Kit	ELISA	142,000	Build-to-order manufacturing
K36	Anti Type IV Collagen NC1 domain Antibodies Test Kit	ELISA	74,000	For human, rat, and monkey. ✖ Build-to-order manufacturing
K61	Collagen Stain Kit	80 test	28,000	This can be quantified simultaneously stained collagen and non-collagen protein. Use for fibrosis research such as a blood vessel wall, liver and kidney.

K62	Collagen stain Red	10ml × 3	24,000	Use for staining of pathological tissue sections.
K71	Human Anti Type II Collagen Antibodies Test Kit	ELISA	64,000	Accompanying antigen / human type II collagen (II C) This can be used to the characteristics analysis of autoantibodies by use with II C of various speacies (K45-49).
K72	Rat Anti Type II Collagen Antibodies Test Kit	ELISA	57,000	Accompanying antigen / bovine II C This can be used to the characteristics analysis of autoantibodies by use with II C of various origin
K73	Mouse Anti Type II Collagen Antibodies Test Kit	ELISA	57,000	Accompanying antigen / bovine II C This can be used to the characteristics analysis of autoantibodies by use with II C of various origin
K77	Type II collagen detection kit	ELISA		Targeting culture supernatant, ELISA.
K79	Type 4 collagen detection kit	ELISA	74,000	For research on diabetic nephropathy and pathological condition mode

FITC

No.	Name	Quantity	MRP	Outline
K21	Type I Collagen FITC (Liquid), Bovine	10ml	32,000	FITC-labeled bovine type I collagen substrate
K21MATU	Type I Collagen FITC (Lyophilized powder), Bovine	10mg	32,000	FITC-labeled bovine type I collagen substrate
K22	Type II Collagen FITC (Liquid), Bovine	5ml	24,000	FITC-labeled bovine type II collagen substrate